HEATTM Collection Smart Pliers

By ThermalOrth Philosophy System

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HEATTM Collection Smart Pliers

By ThermalOrth Philosophy System

The **Heat Collection Smart Pliers®** are the most technologically advanced expression of instruments designed to modify aligners.

With the **Heat Collection** pliers you can create buttons for inter-maxillary elastic tractions, anterior and posterior bites, to make space inside the aligner to speed up dental movements, to restore the surface of the aligner by cancelling activation points, buttons, etc. ...

The orthodontist has the most sophisticated instruments to **instantly** modify aligners, thus speeding up treatment times.




For best performance, **Heat Collection** instruments must be thermoregulated with *iThermal*.

Thermoregulate the tip with the flame symbol, only (or in some models the tip with the writing "heat only here").

It is important to slowly modify the plastic and wait 10/15 seconds before removing the plier, by cooling down with air or water. Carefully handle heated pliers to avoid severe burns from improper use.

The patented graduated scale the **Heat Collections** pliers are provided with, is a precise reference for the operator who can modify the opening of the tips and consequently the size of the thermoforming. To make Buttons, we recommend to deactivate the graduated scale by inserting the key into the Allen screw and turning it counterclockwise, until the screw disappears into its seat.

Take a look at the ThermalOrth Philosophy text book code TP-001-401: it is a complete guide with more detailed instructions to properly use the graduated scale and the **Heat Collection, Hybrid Collection and Cold Collection Smart Pliers**.

| | |
|--|---|
| <p>Button Forming step 1</p>  <p>TP-001-071/SM</p> <p>To make the cylinder.</p> | <p>Button Forming step 2</p>  <p>TP-001-072/SM</p> <p>To make the button.</p> |
| <p>Flattening</p>  <p>TP-001-075/SM</p> <p>To cancel activation points and lines, buttons, etc.</p> | <p>Spacing</p>  <p>TP-001-078/SM</p> <p>To make space inside the aligner</p> |
| <p>Bite Plane</p>  <p>TP-001-077/SM</p> <p>To make bite planes in the front.</p> | <p>Rear Bite</p>  <p>TP-001-081/SM</p> <p>To make bite planes at the rear.</p> |
| <p>Torque Large</p>  <p>TP-001-083/SM</p> <p>To control the torque for central incisors and upper laterals of larger size</p> | <p>Maintenance Instructions - The Heat Collection Smart pliers are manufactured by using the Japanese medical stainless steel to guarantee a longer duration over time. They must be washed and dried thoroughly before being sterilized in autoclave at 121°C/250°F. Subsequently they must be lubricated with synthetic oil code 001-115 (only one drop). The main causes of spots and discolorations are due to a lubricant excess and to an incorrect cleaning (washing and drying not properly carried out) before sterilization. Any yellow or brown spots can be removed by rubbing the pliers with extra fine steel wool.</p> <p>Warranty - The pliers are guaranteed for 5 years. The warranty is void in case of improper use, tampering or use of non-original spare parts. The warranty does not cover labor and travel expenses, which will always and in any case be borne by the buyer.</p> |